

VO2 Experiment



Aug 20th ~25th. 2010

Jong Woo Kim

Contents


Shpyrko Research Group

- Introduction
- Background
- Experiment
- Results
- Safety
- Entrance process
- Appendix

Introduction

Shpyrko Research Group

- Objective
 - X-ray induced Metal-Insulator Transition in VO₂ Film
- Members
 - Sebastian, Moses, Jong Woo, Oleg Shpyrko (PI)
- Facilities



S/W : EPICS




Lakeshore 331 temperature controller

Introduction

Shpyrko Research Group


- Facilities

Sample



Fluorescence detector

Attenuation



Helium or vacuum in & out


Introduction

Shpyrko Research Group

- Beamline : 2-ID-D



2-ID-D



Background


Shpyrko Research Group

- Metal-Insulator Transition(MIT)
- T_{MIT} can be modified by the followings
 - X-ray
 - Voltage
 - Stress, doping, electric field, magnetic field

Background

Shpyrko Research Group

- Voltage induced MIT (Kim et al.)




Kim et al., Appl. Phys. Lett., 90, 23515 (2007)

Background

Shpyrko Research Group

- Voltage induced MIT (Kim et al.)




Kim et al., Appl. Phys. Lett., 90, 23515 (2007)


Background

Shpyrko Research Group

- X-ray induced MIT (Previous experiment, Kiryukhin et al.)
 - Photo-induced & field-induced metallic phases are identical
 - X-ray induced MIT is caused by X-ray photoelectrons and secondary electrons generated in collisions, Not photo-induced oxygen diffusion neither oxygen-diffusion mechanism



Oleg group, previous experiment(2009)



Kiryukhin et al. Nature, 386, 813 (1997)

Experiment


Shpyrko Research Group

Day		Action	Purpose	Result	Remark
1 st	1	Scan Au Fluor.	Locate Au contacts		
2 nd	2	T ↗ (RT→330K)			
	3	Measure R w/XR @332.5K	Choose device	Device 4	
	4	T ↗ ↘ (332.5→352.6K→319.9K)			
	5	V ↗ ↘ (0→3→0V) @310K	V-induced	Metalic until V off	
3 rd	6	Apply 3V @RT	V-induced	$\Delta R = 7k\Omega$	
	7	V ↘ ↗ (3→2→3) @ RT	Joule-heating	R non-recovered	
	8	XR on @3V, RT	XR-V-induced	Went metalic	
	9	T ↗ (RT→355K) @ 0, 1 and 3V	V-induced	V↑ →Transition↓	
	10	V ↘ ↗ (3→0→3)w/ XR off @RT	Joule-heating	R non-recovered	
	11	XR on/off @3V, RT	V-induced	Trans. Found (17kΩ→115Ω)	
4 th	12	XR on w/ CCD @RT ?V	Locate Au contacts		
	13	T ↗ ↘ (RT→400→RT)	R back on	R returned (to 30kΩ)	
	14	XR on @RT, 0V	XR-induced	R reduced, (to 28kΩ)	
	15	XR on @RT, 2V	XR-induced under V	R reduced, (to 27.8kΩ)	
	16	XR on/off @RT, 2V	XR-induced under V	R reduced slowly	
	17	V ↗ ↘ @RT			
5 th	18	T ↗ ↘ (RT→400→RT)	R back on	Return to 30kΩ	
	19	XR on/off @RT, 3V	XR-induced under V		
	20	T ↗ "tail" area w/ XRD	Structural Phase Trans.	To be analyzed	
	21	T ↘ ↗ (RT 345K) V ↗ on D1	Max. allowable V	18V	
	22	T ↗ ↘ (345→400→RT)	R back on	R still low	
	23	Apply 1V, XR on/off	XR-induced trans.		
	24	T ↗ (RT→ 400K)@4V	V-induced	V↑ →Transition↓	
	25	T ↗ (RT→400K)& V ↗ (0→5.2V)			

Results

Shpyrko Research Group

- Voltage effect
 - T_{MIT} decreases with the increasing V (#9, 25)



- Voltage-induced transition occurs at room temperature :
@300K, R reduced by $7k\Omega$ @3V, but less than $2k\Omega$ @0V or 2V (#6,14,15)
- X-Ray effect
 - T_{MIT} decreases with X-ray on, and faster when with V (#8,14,15,16,17,23)

Safety

Shpyrko Research Group


- It is impossible to emphasize the importance of safety too much
 - Don't rush, move slowly



Entrance process

Shpyrko Research Group

Shpyrko Research Group



Orange : required everytime

Green : for the first visitor

* : ID required

Reference

Shpyrko Research Group


- Kiryukhin et al. Nature, **386, 813** (1997)
- H.T Kim et al Appl. Phys. Lett., **90, 23515** (2007)

Appendix

Appendix : Experiment Details

Shpyrko Research Group

1st Day



[Sample]

- To attach Thermocouple detector on sample
- To identify cables/pins


[Trouble Shooting]

- Thermocouple : K-type detector → Move outside
- Multi-meter : Wire short found → Split & Soldering

Appendix : Experiment Details

Shpyrko Research Group


2nd Day



Appendix : Experiment Details

Shpyrko Research Group


3rd Day



Appendix : Experiment Details

Shpyrko Research Group


4th Day



Appendix : Experiment Details

Shpyrko Research Group

5th Day



Appendix : Experiment Data


Shpyrko Research Group

[Beam Information]

E=10.1KeV ID Gap=10.16keV ZP : 86.24mm

Title	Value
IC US(XRF)	280k/sec @200nA/V
IC US(XRD)	107k/sec @100nA/V
IC DS(XRD)	128k/sec @500nA/V

[Scan #1] XRD



Appendix : Experiment Data

Shpyrko Research Group

[Resistance]


Connection	Resistance	Connection	Resistance
1A-1B	35.6k Ω	1A-2A	35.5k Ω
2A-2B	56.2k Ω	2A-3A	43.7k Ω
3A-3B	70.6k Ω	3A-4A	36.0k Ω
4A-4B	43.5k Ω	4A-5A	41.5k Ω
5A-5B	51.3k Ω		

Connection	Resistance	Connection	Resistance	Connection	Resistance
1B-2B	54.8k Ω	1A-2B	89.9k Ω	3A-4B	63.6k Ω
2B-3B	46.5k Ω	2A-1B	97.0k Ω	4A-3B	59.4k Ω
3B-4B	30.0k Ω	2A-3B	73.5k Ω	4A-5B	90.9k Ω
4B-5B	43.3k Ω	3A-2B	72.9k Ω	5A-4B	43.7k Ω


Appendix : Experiment Data

Shpyrko Research Group


[Scan #24] Au Fluor.



[Scan #25] V Fluor.



[Scan #25] Au Fluor.



Zoom in

Both Au & V Fluorescence Seem to reduce
When going +Z (motor move upstream/
Scanning sample down stream)


[Scan #26]

Test using temp controller as dummy motor
zoom in 100um X 100um scan

Appendix : Experiment Data

Shpyrko Research Group


[Scan #27] Au Fluor.



- ← D1
- ← D2
- ← D3
- ← D4
- ← D5

X=6.0 Y=652; Sample slightly rotated

[Scan #27] V Fluor.




R=36.5k Ω , T=300.0K Center on D1 & 30um X30um scan

[Scan #28]

Aborted : Forgot to open shutter

[Scan #29] Au Fluor.



So far, resistance fluctuation
correspond to temperature fluctuation

Appendix : Experiment Data

Shpyrko Research Group

[Scan #30 & 31]


Test CCD trigger & Fluor. time is parallel to CCD time

[Scan #32]


Image device #1 Loop device using temp. controller as dummy CCD & Fluor. Expos. time = 20 sec

16um X 16 um, T=300.2K

[Scan #32] Au Fluor.



[Scan #33] Au Fluor.



[Scan #34]

Changing Temp to 320K. Not collecting CCD data
- Aborted (Forgot to change line & preset time from 20sec)

[Scan #35]

40um X 40um. Live time : 1sec
First two lines of scan are blank because the shutter was closed
Bumped to T up to 330K

[Scan #36]

100um X 100um. T~330K. Drift

Appendix : Experiment Data

Shpyrko Research Group

[Scan #36]

100um X 100um. T~330K


Drift – sample x went from 12 to 9.42

Sample z from 632 to 582


[Scan #37] Au Fluor.

40um X 40um. 41 points. T=330.2K


R=12.59 kΩ



[Scan #44] Au Fluor.



[Scan #44] Resistance



[Scan #38~40] Test

[Scan #41] 1D scan on right contacts to find top one

[Scan #42] 1D scan along top contact (Z=651, -20<X<30)

[Scan #43] Do #42 again

[Scan #44] 2D scan centered in top gap (X=15, Z=651)
20umX20um

Large change does not correlate with temperature fluctuations

Check sample correction R=12.9kΩ


Temp 332.5K always staying below when approaching

Appendix : Experiment Data

Shpyrko Research Group

[Scan #45] Aborted

[Scan #46] Same as #44 with Flour. Real time set 50
R=11.6k Ω



Resistance correlated to temperature
No correlated to X-ray exposure b/w contacts

[Resistance of device @T=332.5]

Device	Resistance
D1	11.54k Ω
D2	13.55k Ω
D3	20.45k Ω
D4	13.00k Ω

[Scan #47, #48]

Dummy scan to record temp. & resistance

Device	Resistance
D1	11.55k Ω
D2	13.47k Ω
D3	20.45k Ω
D4	12.24k Ω

[Scan #49]

Dummy scan 1001 points to record temp & resistance while illuminating sample without OSA

[Scan #50, #51] Realign OSA X, Y

[Scan #52, #53]

Sample Z along right electrodes to find 4th device

Appendix : Experiment Data

Shpyrko Research Group

[Scan #61]

2D scan over device 4

R=12.27k Ω T=332.5K. Beam jumped 4/5 though

[Scan #63]

Heat through transition 332.0K \rightarrow 342.6K

[Scan #64]

Same as #63 w/ 342.6K \rightarrow 352.6K

[Scan #65]

Temp. sweep 320K \rightarrow 355K. 2°/min outside transition, down to 0.2 °/min

Use OSA Y to collect